

## STK Spray validation guide, from STK Sperm Tracker™ range according to the COFRAC recommendations

Preamble: This technology has already been validated by COFRAC in several French laboratories.

### Methodology:

- **Presentation of the technique, apparatus and procedure, field of application and purpose of the verification/method validation**
- **Determination of performance criteria (parameters) to be checked:**
  - Sensitivity, repeatability, intermediate fidelity
  - Specificity and interactions with other fluids
  - Performance on different supports
- **Determination of acceptable specifications or limits (objectives to be achieved) of these criteria**
- **Bibliographic verification**
- **Experimental design and experimental implementation in the laboratory**
- **Compilation and statistical processing of the data obtained**
- **Conclusion and decision on the operational validation of the technique, with regard to the specifications (acceptable limits) initially set**

### Design of experiment and experimental implementation in the laboratory

1. **Sensitivity, repeatability and intermediate fidelity test: Sample preparation: 50 µL deposits**
  - a. Pure semen
  - b. Diluted semen (1/2, 1/5, 1/10, 1/100, 1/1000)
  - c. Negative control: water

➔ Deposit in **triplicate** of **pure sperm** and **different dilutions of sperm** on the same support such as glass benchtop or KIMTECH type paper, then drying for at least 4 hours.

#### 1<sup>st</sup> deposit

NC	Pure semen	Semen 1/2	Semen 1/5	Semen 1/10	Semen 1/100	Semen 1/1000
----	------------	-----------	-----------	------------	-------------	--------------


#### 2<sup>nd</sup> deposit

NC	Pure semen	Semen 1/2	Semen 1/5	Semen 1/10	Semen 1/100	Semen 1/1000
----	------------	-----------	-----------	------------	-------------	--------------


#### 3<sup>rd</sup> deposit

NC	Pure semen	Semen 1/2	Semen 1/5	Semen 1/10	Semen 1/100	Semen 1/1000
----	------------	-----------	-----------	------------	-------------	--------------

➔ Deposit in **triplicate** of **pure sperm** and **diluted to 1/10** on the same support such as glass benchtop or KIMTECH type paper, then drying for at least 4 hours.



➔ Deposit in **triplicate** of **pure sperm** and **diluted to 1/10** from 3 different batches of sperm on the same support such as the glass bench or kimtech type paper, then drying for at least 4 hours. The experiment will last for 3 days knowing that each day, a different operator will spray the solution on the three deposits (negative control, pure semen and semen diluted to 1/10) for each batch (9 deposits in total per day for one operator).



## 2. Specificity and interactions test: Sample preparation: 50 µL deposits

- a. Pure semen
- b. Pure saliva
- c. Pure urine
- d. Pure blood
- e. Pure faeces
- f. Pure semen and pure blood
- g. Pure sperm and pure saliva
- h. Pure semen and pure urine
- i. Pure semen and pure faeces
- j. Negative control

➔ Deposit in **triplicate** of the **different fluids** on the same support such as the glass bench or paper type KIMTECH, then drying for at least 4 hours

### 1<sup>st</sup> deposit

NC	Pure semen	Pure blood	Pure saliva	Pure urine	Pure faeces	Pure semen + pure blood	Pure semen + pure saliva	Pure semen + pure urine	Pure semen + pure faeces
----	------------	------------	-------------	------------	-------------	-------------------------	--------------------------	-------------------------	--------------------------

### 2<sup>nd</sup> deposit

NC	Pure semen	Pure blood	Pure saliva	Pure urine	Pure faeces	Pure semen + pure blood	Pure semen + pure saliva	Pure semen + pure urine	Pure semen + pure faeces
----	------------	------------	-------------	------------	-------------	-------------------------	--------------------------	-------------------------	--------------------------

### 3<sup>rd</sup> deposit

NC	Pure semen	Pure blood	Pure saliva	Pure urine	Pure faeces	Pure semen + pure blood	Pure semen + pure saliva	Pure semen + pure urine	Pure semen + pure faeces
----	------------	------------	-------------	------------	-------------	-------------------------	--------------------------	-------------------------	--------------------------

### 3. Performance test: Preparation of the support to be tested

- a. Smooth (sink, toilet)
- b. Rough (concrete, wall)
- c. Car interior
- d. Waterproof (glass bench)
- e. Natural (wood)
- f. Plastic (hairbrush, pen etc.)

➔ **Triplicata** deposit of **pure semen** and the **different dilutions** of semen (1/10 and 1/100) on each of these supports, then drying for at least 4 hours.

	1 <sup>st</sup> deposit				2 <sup>nd</sup> deposit				3 <sup>rd</sup> deposit			
Smooth/hard material = bathroom sink	NC	Pure semen	Semen 1/10	Semen 1/100	NC	Pure semen	Semen 1/10	Semen 1/100	NC	Pure semen	Semen 1/10	Semen 1/100
Rough material = concrete, roughcast wall	NC	Pure semen	Semen 1/10	Semen 1/100	NC	Pure semen	Semen 1/10	Semen 1/100	NC	Pure semen	Semen 1/10	Semen 1/100
Car interior	NC	Pure semen	Semen 1/10	Semen 1/100	NC	Pure semen	Semen 1/10	Semen 1/100	NC	Pure semen	Semen 1/10	Semen 1/100
Waterproof material = glass bench	NC	Pure semen	Semen 1/10	Semen 1/100	NC	Pure semen	Semen 1/10	Semen 1/100	NC	Pure semen	Semen 1/10	Semen 1/100
Natural material = wood, leaf	NC	Pure semen	Semen 1/10	Semen 1/100	NC	Pure semen	Semen 1/10	Semen 1/100	NC	Pure semen	Semen 1/10	Semen 1/100
Plastic material = hairbrush	NC	Pure semen	Semen 1/10	Semen 1/100	NC	Pure semen	Semen 1/10	Semen 1/100	NC	Pure semen	Semen 1/10	Semen 1/100

## 4. Preparation of the STK Spray solution

- a. Open a STK Spray pouch and pour it into a sprayer
  - The sprayer must allow to vaporize in fine misting
- b. Mix the powder with 100 mL of demineralized water
- c. Homogenize the solution until the powder dissolves completely (about 45 seconds)
- d. Prime the sprayer by pressing the trigger 2-3 times towards a trash can or sink

## 5. Use of different UV lamps if possible to compare the difference in lamp performance with respect to signal detection.

**Contact AXO Science to find out if your lamp model is suitable.**

**365 nm UV lamp used and recommended by AXO Science: VILBER VL-6.L.**

- a. 365 nm UV lamp → Characteristics to be respected:
  - Portable or plugged into the power supply( at your convenience)
  - Specificity with sperm once STK Spray sprayed
  - Signal to background noise ratio (very little interference): specific and intensity
  - Signal expression (the stain) + easy to distinguish compared to areas devoid of sperm
  - High power of the UV lamp
  - If neon tube lamp: length of the neon (neon tube quite long with a larger surface to be observed)
- b. Protect your eyes well with transparent UV-resistant glasses

## 6. Tests realizations

- a. Being in the dark
- b. Spray the different supports using the sprayer containing the solution by keeping it vertical
- c. After waiting between 45 seconds to 1 minute, illuminate the samples with the 365nm UV lamp (or UV lamps if lamp comparison study):
  - for sensitivity test
  - for specificity test
  - for performance test
- d. Take pictures of the results obtained with a camera whose settings are fixed for each element to be photographed.

***Example for a CANON EOS 2000D camera that is attached to a steel box allowing you to be as much as possible in the dark:***

- ***Focal: F5.6***
- ***ISO: 3200***
- ***Shutter speed: 1/20th of a second***

## 7. Results compilation and data processing

## 8. Conclusion and decision on the operational validation of the technique, against the specifications (acceptable limits) initially set.

